

STAPLE HEIGHT COMPARISON

Select a stapler reload with the appropriate staple size for the tissue thickness

ENDO GIA™ with TRI-STAPLE™

AEON™ ENDOSTAPLER

Gray
2.0mm / 0.75mm

Gray
2.0mm / 0.75mm

Tan
2.0-3.0mm / 0.75-1.25mm

White
2.5mm / 1.0mm

Purple
3.0-4.0mm / 1.25-1.75mm

Orange (for a drier staple line)
3.25mm / 1.5mm

Purple (for more perfusion)
4.0mm / 1.8mm

Black
4.0-5.0mm / 1.75-2.25mm

Black
5.0mm / 2.2mm

STAPLE HEIGHT COMPARISON

Select a stapler reload with the appropriate staple size for the tissue thickness

ECHELON FLEX™ STAPLER

AEON™ ENDOSTAPLER

Gray
2.0mm / 0.75mm

Gray
2.0mm / 0.75mm

White
2.5mm / 1.0mm

White
2.5mm / 1.0mm

Blue
3.5mm / 1.5mm

Orange
3.25mm / 1.5mm

Gold
3.8mm / 1.8mm

Purple
4.0mm / 1.8mm

Green
4.1mm / 2.0mm

Black
4.2mm / 2.3mm

Black*
5.0mm / 2.2mm

*The AEON Black Reload uses a 15mm trocar