


STAPLE HEIGHT CHART



Echelon Flex™ Staple Height (mm)		
	OPEN	CLOSED
Gray	2.0	0.75
White	2.5	1.0
Blue	3.5	1.5
Gold	3.8	1.8
Green	4.1	2.0
Black	4.2	2.3

Endo GIA™ Staple Height (mm)		
	OPEN	CLOSED
Gray	2.0	0.75
White	2.5	1.0
Blue	3.5	1.5
Green	4.8	2.0

Endo GIA™ with Tri-Staple™ Staple Height (mm)		
	OPEN	CLOSED
Gray	2.0*	0.75
Tan	2.0, 2.5, 3.0	0.75, 1.0, 1.25
Purple	3.0, 3.5, 4.0	1.25, 1.5, 1.75
Black	4.0, 4.5, 5.0	1.75, 2.0, 2.25

*Per Endo GIA™ brochure, Extra-Thin/Vascular reload contains three rows of 2.0mm staples (non-varied height) and contains all of the other features and benefits of an Endo GIA™ Reload with Tri-Staple™ Technology.