

Woop

WOUND PROTECTOR



Maximum Protection Minimum Infection



ankalaps®
www.ankasurgical.com

ISO 13485:2016

CE
2195



White ring is easy to roll down to adjust to variable thicknesses of the abdominal Wall

Strong Polyurethane material offers optimal wound protection

Flexible distal blue ring is inserted through the incision and easily removed at the end of the procedure

GENERAL

Inguinal Hernia Repair (XS, S), Thyroidectomy (XS, S) Appendectomy (S, M), Splenectomy (L, XL)

BARIATRIC

Lap Gastric Bypass (XS, S), Open Gastric Bypass (L, XL)

COLON & RECTAL

Lap Colectomy (S, M Laparoscopic System)
Open Colectomy (L, XL)

CARDIOTHORACIC

Video-Assisted Thoracoscopic Surgery (VATS) (XS, S) Mitral Valve Repair/Replacement (S, M)
Thoracotomy (S, M)

ORTHOPAEDIC

Total Shoulder Arthroplasty (S), Total Hip Arthroplasty (S, M, L)

OB/GYN

Postpartum Tubal Ligation (XS), Bilateral Salpingo Oophorectomy (XS, S), Mini-Laparotomy (S, M) Myomectomy (S, M), Total Abdominal Hysterectomy (S, M, L), Cesarean Section (L, XL)

BREAST

Lumpectomy (XS, S), Mastectomy (S, M), Sentinel Lymph Node Biopsy (XS, S)

SQUARE MODEL



INCISION PROTECTOR RETRACTOR					
	INCISION LENGTH	HARD / ROUND	HARD / SQUARE	SOFT / ROUND	SOFT / SQUARE
SIZE		PRODUCT CODE			
XXS	1.0 - 3.0 cm	AIPXXS-HR	AIPXXS-HS	AIPXXS-SR	AIPXXS-SS
XS	2.0 - 4.0 cm	AIPXS-HR	AIPXS-HS	AIPXS-SR	AIPXS-SS
S	2.5 - 6.0 cm	AIPS-HR	AIPS-HS	AIPS-SR	AIPS-SS
M	5.0 - 9.0 cm	AIPM-HR	AIPM-HS	AIPM-SR	AIPM-SS
L	9.0 - 14.0 cm	AIPL-HR	AIPL-HS	AIPL-SR	AIPL-SS
XL	11.0 - 17.0 cm	AIPXL-HR	AIPXL-HS	AIPXL-SR	AIPXL-SS
XXL	14.0 - 20.0 cm	AIPXXL-HR	AIPXXL-HS	AIPXXL-SR	AIPXXL-SS

INCISION PROTECTOR RETRACTOR WITH LAPAROSCOPIC COVER				
	INCISION LENGTH	HARD / ROUND		HARD / SQUARE
SIZE		PRODUCT CODE		
XS	2.0 - 4.0 cm	4 PORT	AIPXS-HRC	AIPXS-HSC
S	2.5 - 6.0 cm	1 PORT	AIPS-HRC	AIPS-HSC
M	5.0 - 9.0 cm	1 PORT	AIPM-HRC	AIPM-HSC